







Photo	Item No	Product Names	Specification
	YSX200VET Plus	Digital 16KW Veterinary X-Ray Machine	<b>Toshiba tube</b> <b>16KW/200mA</b> , suitable for all parts of large and small animals.
	YSX200vet	16KW Veterinary X-Ray Machine	<b>Toshiba tube</b> <b>16KW</b> , suitable for all parts of large and small animals.
	YSX050-A	Portable veterinary x-ray machine	5KW/100mA touch screen with mobile stand
	YSX050-A	Portable veterinary x-ray machine	5KW/100mA button with mobile stand
	YSX050-B	Portable veterinary x-ray machine with table	5kw/100mA touch screen with table
	YSX050-C	<b>Digital</b> Portable veterinary x-ray machine	5KW/100mA portable x-ray machine mobile stand button control panel flat panel detector computer veterinary software

Photo	Item No	Product Names	Specification
	YSX1509	Lead apron	0.35mmPB
	YSX1513	Lead apron	0.35mmPB
	YSX1520	X-ray lead gloves	0.35mmPB
	YSX1603	Lead Glasses	0.50mmpb
	YSX1536	Lead sheet for x-ray room protection	<b>Before we quote, we need to know your xray room size and thickness of lead sheet you require.</b>
	YSX1525	Leaded door for x-ray room	<b>Please let us know your xray room door hole (without frame) size, then we can quote the exact price.</b>
	YSX1613	Leaded glass	<b>For vet x-ray room, normally a glass sized 40CM*40CM is enough. If you require other size, please let us know.</b>
	YSB-W1	Wireless veterinary convex probe	<b>With a HUAWEI mobile phone</b>
	YSB-W2	Wireless veterinary linear probe	<b>With a HUAWEI mobile phone</b>
	YSB180V <b>HOT!</b>	Portable veterinary ultrasound scanner	5.0MHz micro-convex probe
	YSB5600V <b>HOT!</b>	Portable veterinary ultrasound scanner	5.0MHz micro-convex probe
	YSB-K10V	Portable veterinary ultrasound scanner	6.5MHz rectal linear probe With WIFI function
	YSB5000V	Portable veterinary ultrasound scanner	6.5MHz rectal linear probe for large animals










Photo	Item No	Product Names	Specification
	YSX050-C	<b>Digital</b> Portable veterinary x-ray machine	5KW/100mA portable x-ray machine mobile stand <b>touch screen</b> <b>flat panel detector</b> computer veterinary software
	YSX50M	50mA mobile vet mobile x-ray machine	50mA, for arms and legs of large animal, and for all parts of small animal
	YSX100VET	100mA mobile veterinary x-ray unit	<b>100mA</b> <b>Table is optional</b>
	YSX100VET	100mA Digital mobile veterinary x-ray unit	<b>5kW / 100mA Digital</b> <b>Table is optional</b>
	YSX1501	Automatic film processor	Developing Capacity: 80p/h (14"*17"); 145p/h(10"*12") video link: <a href="https://www.youtube.com/watch?v=LgHcMoRU9IY">https://www.youtube.com/watch?v=LgHcMoRU9IY</a>
	YSX1702	Dark roomlight	Single color
	YSX1707	Cassette	8"*10"
			11"*14"
			12"*15"
	YSX1615	Intensifying screen (green)	8"*10"
			11"*14"
			12"*15"
	YSX1619	x-ray film (green)	8"*10"*100 sheets
			11"*14"*100 sheets
			12"*15"*100 sheets

Photo	Item No	Product Names	Specification
	YSB5100V	Palm type veterinary ultrasound scanner	6.5MHz rectal linear probe for large animals
	YSB5200V	Portable veterinary ultrasound scanner	6.5MHz rectal linear probe for large animals
	YSB-U3V	Portable veterinary ultrasound scanner	With sector probe
	SonoScape E1V	Portable veterinary B/W ultrasound scanner	With micro-convex probe
			With rectal & linear probes
	SonoScape S2V <b>HOT!</b>	Veterinary portable color doppler ultrasound	With one micro-convex and one rectal probe
	YSB-V30 VET	Veterinary portable color doppler ultrasound	With one micro-convex and one rectal probe
	SonoScape S11V	Trolley color doppler ultrasound	With one micro-convex and one rectal probe
	YSPM400V <b>HOT!</b>	Veterinary patient monitor	8" screen ECG /RESP /SpO2 /NIBP /TEMP /Lithium Battery
	YSPM500V	Veterinary patient monitor	12.1" large-sized screen, anti-dazzle with 4° incline ECG /RESP /SpO2 /NIBP /TEMP /Lithium Battery

	YSPM200V	Veterinary vital sign monitor	7"high-precision color LCD display ECG /RESP /SpO2 /NIBP /TEMP /Lithium Battery
	YSPO100V	Veterinary Handheld Pulse Oximeter	SpO2 /Pulse Rate /Hanging Strap
	YSECG300V	Veterinary 3-channel ECG machine	1. 3 channel ECG machine with interpretation 2. 5.0" color TFT LCD display 3. Simultaneous 12 leads acquisition and 1, 1+1, 3 channel (Manual/Auto) recording with high resolution thermal printer 4. Manual/Auto working modes
	YSECG-i3V	Veterinary 3-channel ECG machine	5.6" high resolution LED screen
	YSBP80V	Veterinary blood pressure monitor	High-definition color LCD display, strong visibility. Data storage function, up to 100 groups of data can be stored, by which measured data can be reviewed.
	YSDBP100V	Veterinary doppler ultrasound blood pressure system	1. Measuring range: 0mGs~199.9mGs 2. Frequency of the probe: 8MHz 3. Specification of the blood pressure cuff Range: No.1 4.2cm~7.1cm No.2 5.4cm~9.1cm No.3 6.9cm~11.7cm No.4 9cm~13cm No.5 12cm~16cm 4. Rechargeable battery voltage:7.4V ,1500mAh 5. AC input: 100~240V 50/60Hz DC output
	YSDBP200V	Veterinary doppler ultrasound blood pressure system	1. Measuring range:0mGs~199.9mGs 2. Frequency of the probe: 8MHz 3. Specification of the blood pressure cuff Range: No.1 4.2cm-7.1cm No.2 5.4cm-9.1cm No.3 6.9cm-11.7cm No.4 9.0cm-13.0cm No.5 12.0cm-16.cm No.6 17.0cm-22.0cm 4. Rechargeable battery voltage: 7.4V, 1500mAh 5. AC Input: 100-240V 50/60Hz DC Output
	YSTE880V	Animal blood analyzer	Fully automatic Single channel 35 tests/h

Photo	Item No	Product Names	Specification
	YSTE320V <b>HOT!</b>	Veterinary hematology analyzer	Double channels 60 tests/h <b>10.4 inch touch screen</b>
	YSTE-DF50V <b>TOP!</b>	Veterinary hematology analyzer	Double channels 60 tests/h <b>5-part</b>
	YSTE-R5V	Reagent for YSTE-DF50V	DIL-C for VET (20L) LYC-1 for VET (200ml) LYC-2 for VET (500ml) CLE-P (50ml)
	YSTE-21BV <b>HOT!</b>	Semi-auto veterinary chemistry analyzer	5-inch touch screen and keypad Wavelength: 340、380、405、505、546、578、620nm
	YSTE168V	Semi-auto veterinary chemistry analyzer	Wavelength: 7 standard filters, 340, 405, 505, 546, 578, 620 670nm, plus one more optional filter
	YSOT-LED50MD	LED Operation lamp	LED, Illumination Intensity (lx): 30000-160000, adjustable <b>With battery</b>
	YSVET8006	Veterinary OT table	201 stainless steel L100cm*W60cm*H45--105cm
	YSVET0506	V-shap vet operation table	304 Stainless steel <b>(with heater and grid filter)</b> Size:1300*600*450/950 mm

Photo	Item No	Product Names	Specification
	YSVET0506	V-shap vet operation table	304 Stainless steel <b>(without heater and grid filter)</b> Size:1300*600*450/950 mm
	YSVET0506	V-shap vet operation table	304 Stainless steel <b>(with grid filter)</b> Size:1300*600*450/950 mm
	YSVET0506	V-shap vet operation table	304 Stainless steel <b>(with heater)</b> Size:1300*600*450/950 mm
	YSVET0506	Flat-top vet operation table	304 Stainless steel <b>flat top</b> Size:1300*600*450/950 mm
	YSVET0507	Veterinary surgery table	High quality 304 stainless steel Size: 1500*600*760/1050mm
	YSVET0504	Animal <b>electric</b> operating table	1、 Table height: 760—1060mm 2、 Table dimension: 1400*650mm 3、 Table gradient : ±15° 4、 Adjustable angle of table board: 0°-45° 5、 The range of temperature: 0°—50° 6、 Power: 220V
	YSVET3101	Veterinary dissectign table	High quality stainless steel Size: 1400mm*600mm *830mm
	YSVET1600	Veterinary dissectign table	High quality 304 stainless steel Size: large: 1600*650*800mm Small: 110*550*800mm
	YSVET3102	Veterinary dissectign table for small animals	High quality 304 stainless steel Size: 690mm*300mm*100mm

Photo	Item No	Product Names	Specification
	YSVET1360	Veterinary stretcher cart	High quality 304 stainless steel 1360mm*600mm*850mm
	YSVET-DJ02 <b>NEW!</b>	Veterinary stretcher	120*55CM
	YSVET5102	Surgical instrument table	High quality 304 stainless steel 800mm*500mm*830mm
	YSVET-DZ01	Animal Injection Cage	S: 700*360*500mm M: 970*420*680mm L: 1200*460*780mm
	YSVET-BDJ	Pet restraint table	Quality 304 stainless steel 600mm*450mm*800mm 
	YSVET-SGQ01	Leash tying hook	Base: 9*9cm
	YSVET-ICU <b>NEW!</b>	<b>Pet ICU Chamber &amp; Infant Incubator</b>	Integrated the functions of pet temperature recovering, infusion, oxygen supply, phototherapy, nebulizing therapy, emergency rescue, hospitalized observation, etc. It can also be an infant pet incubator
	YSVET-UC1801	Pet Incubator	Accurate Temperature Control—Being Equal Everywhere Accurate Humidity Control—Healthy and Balanced Negative-ion Generation—Truly Auxiliary Medical Effect Sterilization Function—Free From Cross-contamination Medical Atomization Treatment Function—Making the Treatment More Convenient Carbon Dioxide Concentration Monitoring System—A Defense Line to Guard Life ICU Illumination Function—Creating A Comfortable Therapeutic Environment



Photo	Item No	Product Names	Specification
	YSVET610A	Animal oxygen supply cage	61*61*70cm
	YSVET610B	Veterinary cage	High quality 304 stainless steel Size: 610*700*610/1410mm
	YSVET0510 <b>TOP!</b>	Combination cage	Stainless steel, size: 2400*700*2140mm, including two big cages, three middle cages, five small cages and one oxygen cage
	YSVET2440	Pet display cage	Stainless steel, size: 2440*700*2180mm, For separate cages at each layer
	YSVET1220	Stainless steel pet boarding cage	# 304 stainless steel External size: 1220*700*1580mm (L*D*H) 1. Upper cage: 610*610*700mm (L*D*H) 2. Lower cage: 610*820*700mm (L*D*H) 3. Waste tray: 690*540*45mm (L*W*H)
	YSVET1220A	Veterinary in-patient cage	# 304 stainless steel External size: 1220*700*2180mm (L*D*H) Top layer: 610*610*700mm (L*D*H) Middle layer: 610*610*700mm (L*D*H) Bottom layer: 610*820*700mm (L*D*H) Waste tray: 690*540*45mm (L*W*H) <b>Infusion pump is optional</b>
	YSVET1220B	Stainless steel pet boarding cage	# 304 stainless steel External size: 1220*700*1570mm (L*D*H) 1. Upper cage: 610*610*700mm (L*D*H) 2. Lower cage: 610*820*700mm (L*D*H) 3. Waste tray: 690*545*10mm (L*W*H)
	YSVET1220C	Stainless steel pet cage <b>with infrared light</b>	# 304 stainless steel External size: 1220*700*1570mm (L*D*H) 1. Upper cage: 610*610*700mm (L*D*H) 2. Lower cage: 610*820*700mm (L*D*H) 3. Waste tray: 690*545*10mm (L*W*H)
	YSVET1000	Pet cage with spring board	With spring board High quality 304 stainless steel Size: 1000*810*650 1. External size: 1000*810*650mm 2. Upper cage: 1000*810*650mm Lower cage: 1000*810*650mm 3. Springboard: 580*250*10mm 4. Left door: 670*750*8mm 5. Right door: 270*750*8mm

Photo	Item No	Product Names	Specification
	YSVET-ZGB01 <b>NEW!</b>	Block-out cloth/shade cloth	<b>To be tailor-made with cages</b>
	YSVET-BD04 <b>NEW!</b>	Pet ice mat	Size: 65*50cm
	YSVET-DRB01 <b>NEW!</b>	Thermostat carbolic-crystal therapy board	Size: 50*50*1.2cm Power: 60W Power supply: 220V Lifespan: 10 years
	YSVET8105	Pet transfer cage	Full 304 stainless steel 1220mm*690mm*960mm
	YSVET0508 <b>HOT!</b>	Stainless vet grooming table	#304 Stainless steel, the length of the tap can be adjustable Machine size: 140*60*88cm
	YSVET0508A	Stainless vet treatment table	#304 Stainless steel Machine size: 100*60*75/105cm
		Pet spa tub	Made of quality acrylic & #304 stainless steel Size: 1100*900*850mm Computer control with bubbling and spa function Pure acrylic, environmental friendly and healthy With resin structure to make it more durable Fibreglass ensure pressure-proof Sleek shape, seamless joint Water draining by hand wheel switch Interior tub is water-tight and water-shedding With pumping faucet

Photo	Item No	Product Names	Specification
	YSVET-CX130 A	Spa tub	Fully 304 stainless steel 1300*740*900/1300mm
	YSVET-CX130 B	Pet cleaning tub with hair dryer	Fully 304 stainless steel 1300*740*900/1300mm
	YSVET-CX130 C	Pet cleaning scrub	Fully 304 stainless steel 1300*740*900/1300mm
	YSVET-CX130 D	Pet cleaning scrub	Fully 304 stainless steel 1300*740*900/1300mm
	YSVET-CX130 E	Pet cleaning tub with stepping board	Fully 304 stainless steel 1300*740*900/1300mm
	YSVET-QX620	Cleaning sink	Fully 304 stainless steel 620*500*800mm
	YSVET-QX620A	Cleaning sink	Fully 304 stainless steel 620*500*800mm
	YSVET-QX620B	Cleaning sink	Size: 62*50*120cm

Photo	Item No	Product Names	Specification
	YSVET-QX1200	Double-person wash sink	Fully 304 stainless steel 1200*600*1700mm
	YSVET10818	Portable veterinary hair drier	Voltage: 220V-240V 50hz Power: 2800W 3.7HP Air speed: 27000fpm-36000fpm Nois: 68DB Heat: 43℃-60℃
	YSVET10901L	Standing pet hair drier	Voltage: 220V-240V 50hz Power: 2800W 3.7HP Air speed: 27000fpm-36000fpm Noise: 68DB Heat: 43℃-60℃
	YSVET8908	Hair clipper	Battery type: Ni-MH/Ni-cd Battery spec.: AA700-1200mah 2.4V Charger input: AC220-230V 50Hz Charger output: C3V300-800mA Maximum power: 3-10W Charging time: 6-8h
	YSVET-MY8002	Veterinary grooming table	L90cm*W60cm*H75cm
	YSVET-MY8003	Veterinary grooming table	L120cm*W60cm*H75cm
	YSVET-MY8005	Veterinary grooming table	L120cm*W60cm*H45--105cm
	YSVET-MY1001	Veterinary grooming table	Table-top diameter: 800 mm Height adjustment: 760 mm to 830 mm

Photo	Item No	Product Names	Specification
	YSVET-MY1002	Veterinary grooming table	1000 mm * 600 mm * 760/830mm
	YSVET-MY829	Veterinary grooming table	<b>LED lighted table top!</b> 120×60×25.5-95
	YSVET-MY829	Veterinary grooming table	<b>LED lighted table top!</b> 120×60×25.5-95
	YSTZC-040	Wighing scale for small animals	50g~20kg
	YSTZC-01	Animal scales	Material: pattern steel Battery: can charge Biggest weighing: 1t - 10 tons, can be customized
	YSVET-TZC150	Animal scales	Full 304 stainless steel Loading capacity: 150KG
	YSVET-ZDA <b>NEW!</b>	Wall-mounted foldable treatment table	Table top: 110*56cm 

Photo	Item No	Product Names	Specification
	YSVET2101	Veterinary examination & treatment table	High quality 304 stainless steel Size: 1200mm*600mm*830mm
	YSVET2102	Veterinary examination & treatment table	With weighing scale Table-board length and width: 1200mm*600mm Table height: 830mm Fence height: 300mm
	YSVET2103	Veterinary examination & treatment table	High quality 304 stainless steel Size: Large: 1100*600*800mm Small: 900*500*800mm
	YSVET2104	Veterinary examination & treatment table	High quality 304 stainless steel Size: 1000*550*800mm
	YSVET2105	Veterinary examination & treatment table	High quality 304 stainless steel Size: 1300*600*450-1000mm
	YSVET2106	Veterinary examination & treatment table	High quality 304 stainless steel Size: Large: 1100*600*800mm Small: 900*500*800mm
	YSVET2107	Veterinary examination & treatment table	High quality 304 stainless steel Size: 920*520*840mm
	YSVET2108	Veterinary examination & treatment table	High quality 304 stainless steel Size: Part one: 1100*520*770mm Part two: 650*300*770mm
	YSVET-BCT01 <b>NEW!</b>	Veterinary table for ultrasound scanner	Size: 110*55*70cm
	YSVET-F06 <b>NEW!</b>	Trolley cart for veterinary equipment	Size: 600*650*1525mm
	YSVET1101	Veterinary infusion table	High quality 304 stainless steel Size: 1200mm*600mm*830mm

Photo	Item No	Product Names	Specification
	YSVET1102	Veterinary infusion table	With heating function Size: 1200mm*600mm*830mm
	YSVET1103	Veterinary infusion table	High quality 304 stainless steel Large: 1100*600*850mm Small: 900*500*850mm
	YSVET1104	Veterinary infusion table	High quality 304 stainless steel Large: 1100*600*780mm Small: 900*550*780mm
	YSVET1105	Veterinary infusion table	High quality 304 stainless steel Large: 1100*600*800mm Small: 900*500*800mm
	YSVET1106	Veterinary infusion table	High quality 304 stainless steel Large: 1100*600*850mm Small: 900*500*850mm
	YSVET1107	Veterinary infusion table	High quality 304 stainless steel Large: 1100*600*780mm Small: 900*550*780mm
	YSVET1108	Veterinary infusion table	High quality 304 stainless steel Size: 1370*600*780mm