

CERES



CERES

Semi-Automatic Coagulometer
Product code 3801100

Principle: Scattered light with percentage analysis
Test Channel: 2
Light source: Durable LED detector
Test time: 20-60 secs
4 sample position
2 reagent position (1 magnetic stirrer included)
Reagent consumption: 25µl for APTT/TT/FIB, 40µl for PT
Memory: 6.000 test results
Interface: RS-232
Display: LCD display
LCD input: Easy operation keypad
Output: internal printer
Weight: 5 kg
Dimensions: 28x31x16 cm

CERES

Coagulomètre Semi-automatique.
Code produit 3801100

Principe: Diffusion de la lumière avec analyse en pourcentage
Canal test: 2
Source de lumière: LED durable
Durée du test: 20-60 secondes
4 position échantillon
2 position réactif
Consommation du réactif: 25µl pour APTT / TT / FIB, 40µl pour PT
Mémoire: 6 000 résultats
Interface: RS-232
Écran: LCD
Input: clavier
Output: imprimante interne
Poids: 5 kg
Dimension: 28x31x16 cm

CERES

Coagulómetro semi-automático.
Referencia producto 3801100

Principio: Dispersión de luz con análisis porcentual
Nº de canales: 2
Fuente de luz: LED de larga duración
Tiempo por test: 20-60 segundos
4 posiciones de muestra
2 posiciones de reactivo
Consumo de reactivo : 25µl for APTT/TT/FIB, 40µl para PT
Memoria: 6.000 test results
Comunicación: RS-232
Pantalla: LCD
Input: Teclado de fácil uso
Output: Impresora interna
Peso: 5 kg
Dimensiones: 28x31x16 cm